



The Beauty of New IGC Applications

Dr. Ralf Dümpelmann, Inolytix AG, Switzerland

ralf.duempelmann@inolytix.com

+41 79 3129327



Exclusive partner of Adscientis for
Germany, Austria and Switzerland



3 Examples ... to inspire




<https://photoartflight.com/2016/09/21/black-swallowtail-butterfly-takeoff/>


Specific surface area (BET) by water



Flat surface investigations – coatings and films



Diffusion coefficients by Zero Length Column (ZLC)



Specific surface area (BET) by water

Question

Interaction with H₂O
(soils, cement)

Idea / concept

H₂O as probe
molecule,
TCD-Detector


Results

IGC-FC: $S_{\text{BET}}(\text{H}_2\text{O})$,
heterogeneity -
AEDF

Hydrophilic and hydrophobic soil
(Univ. of Western Australia)



IGC-FC: saturation of sample with
water-vapour and desorption curve




Specific surface area (BET) by water


Results

IGC-FC: $S_{\text{BET}}(\text{H}_2\text{O})$,
heterogeneity -
AEDF

BET-regression: highly linear!



Calculation of adsorption/desorption isotherm: adsorbed amount vs. p/p_0



Code	Q_0 [$\mu\text{mol/g}$]	S_{BET} [m^2/g]	C_{BET} [n.u.]
WET	$121,4 \pm 2,1$	$8,7 \pm 0,2$	$11,6 \pm 1,2$
REP	$111,5 \pm 2,0$	$8,0 \pm 0,1$	$12,0 \pm 0,4$

Specific surface area (BET) by water


Results

Comparison of H₂O (hydropilic) and CH₂Cl₂ (hydrophobic)

Code	Probe	Q ₀ [μmol/g]	S _{BET} [m ² /g]	C _{BET} [n.u.]
WET	Water	121.4 ± 2.1	8.7 ± 0.2	11.6 ± 1.2
REP	Water	111.5 ± 2.0	8.0 ± 0.1	12.0 ± 0.4
WET2	Dichloromethane	21.6 ± 0.2	3.9 ± 0.1	13.9 ± 0.4
REP2	Dichloromethane	22.8 ± 0.2	4.1 ± 0.1	8.9 ± 0.3


Apparent specific surface area $S_{BET}(H_2O) \gg S_{BET}(CH_2Cl_2)$

BUT: only small differences between wettable + repellent



Other causes of reppellency?

2% organic C...



Flat surface investigations - coatings and films

Question


How to measure flat surfaces? (coated metal sheets, films, polymers)

Idea / concept

use a long gas path across the surface (10x10 cm, 29 cm²)

Results

normal IGC-ID and IGC-FC tests




<https://coleofduy.com/wp-content/uploads/2020/04/Adhesive-Film-Market.png>

Flat surface investigations - coatings and films

Results
normal IGC-ID
and IGC-FC tests

- Adhesive films
 - Characterization of the polymer film
 - Characterization of the substrate
 - Compatibilities – and optimizations
 - OPEN ISSUES...

- Coated or treated metal sheets:
 - BET surface characterization, direct!
 - Surface properties
 - Use of special reagents
 - Application note (soon)




IGC on flat surfaces

Example of a steel plate degreased using two processes:

Sample A : Treatment A

Sample B : Surfactant degreasing

Objective : compare the “cleanliness” of both treatments



$$S_{GEO} = 28.8 \text{ cm}^2$$

		$S_{BET} [\text{m}^2]$	S_{BET}/S_{GEO}
Hexanol	A	<u>0.072 ± 0.016</u>	8.0
	B	0.044 ± 0.002	4.9
Hexylamine	A	<u>0.071 ± 0.014</u>	7.9
	B	0.060 ± 0.005	6.7

Hexylamine and hexanol are probes for “cleanliness”:

$$S_{BET} (A) > S_{BET}(B)$$

=> Treatment A more efficient than B

$$S_{BET} (\text{hexanol}) = S_{BET}(\text{hexylamine}) \text{ for A}$$


=> Treatment A: more homogeneous

$$S_{BET} (\text{hexanol}) < S_{BET}(\text{hexylamine}) \text{ for B}$$

=> Treatment B: a part of the basic sites remain inaccessible

Surface heterogeneity: even on flat surfaces

Adsorption Energy Distribution Function (AEDF)



Two kinds of sites (bimodal AEDF)

Significantly less basic sites on B

Slightly less acidic sites on B

**Treatment A is more efficient
to regain adhesive strength of the steel
(especially visible with hexanol)**

Diffusion coefficients by IGC

Question

How to improve current D_{eff} determination?



Idea / concept

use IGC-FC possibility and desorption curve

Results


D in porous samples by any gaseous probe

Currently: van Deemter – variation of flow rates and calculation of HETP (height equivalent to a theoretical plate) **Time consuming!**



Zero Length Column (ZLC)

Saturation at p/p_0



Diffusion coefficients

Results

D_{eff} in porous samples by any gaseous probe

ZLC method:



Conversion of the desorption front and fitting the asymptote in a $\ln(C/C_0)$ plot to obtain D .

$$\ln\left(\frac{C}{C_0}\right) = \ln\left(\frac{2L}{\beta^2 + L(L-1)}\right) - \beta^2 \frac{D \cdot t}{R^2}$$

References

M. Eic, D.M. Ruthven, Zeolites, 8(1), 40–45, (1988)

D. Ruthven, F. Brandani, Adsorption, 11(1), 31–34, (2005)



Example: Diffusion coefficients


Example: single Alumina beads

B12: $d_{50} = 1.71$ mm

B17: $d_{50} = 1.65$ mm

$$D_{\text{toluene}} < D_{\text{cyclooctane}} < D_{\text{octane}}$$


toluene



$$D_{\text{tol}} = 2.14 \cdot 10^{-5} \text{ cm}^2/\text{s}$$

$$D_{\text{tol}} = 2.07 \cdot 10^{-5} \text{ cm}^2/\text{s}$$


cyclooctane



$$D_{\text{cyc}} = 2.85 \cdot 10^{-5} \text{ cm}^2/\text{s}$$

$$D_{\text{cyc}} = 2.50 \cdot 10^{-5} \text{ cm}^2/\text{s}$$

n-octane



$$D_{\text{oct}} = 3.63 \cdot 10^{-5} \text{ cm}^2/\text{s}$$

$$D_{\text{oct}} = 3.87 \cdot 10^{-5} \text{ cm}^2/\text{s}$$

Beauty of many possibilities

- Flexible use of (almost) any vaporizable probe
 - Custom-oriented tasks: corrosive gases, functional groups (double bonds...)
- Saturation at p/p_0 by constant flow with any probe (vaporizable)
- Specific surface area measurements with any probe
- with 22+ probes: good coverage of 3-dimensional HSP sphere
 - HSP: Hansen Solubility Parameters for formulations
- Detectors for organics (FID) and inorganics (TCD)
- Device for flat surfaces

How would you like
to characterize your samples
with 15+ gas probes?



Take your opportunities ...
... and use your imagination!

Demands, ideas or specific questions?

Don't hesitate to ask!

Dr. Ralf Dümpelmann

ralf.duempelmann@inolytix.com

+41 79 312 9327

